



### Getting used to your Clarity Chat

On average it takes 2-3 weeks of daily usage to become fully accustomed to wearing.



Since everyone adjusts to using sound amplifiers at a different rate, we suggest wearing your Clarity Chat for a few hours on the first day, and gradually increasing your usage over the course of the first week.



When you first begin using your device, keep in mind that your brain may be startled when it begins to receive sounds and signals it's been missing. Car engines, keys dropping, doors shutting, typing on a keyboard, or even your own voice may seem loud at first. With consistent usage over time, your brain will adjust to these environmental noises so they are no longer startling and you can begin to tune out these unwanted sounds.

### Operating your Clarity Chat

#### How to turn on your Clarity Chat

Please see diagram 1 on back for battery installation instructions

#### How to wear your Clarity Chat

Please see diagram 2 on back for placement instructions


#### How to use and adjust the sound of your Clarity Chat

##### Programs:

- Program 1** Recommended program for all listening situations
- Program 2** Recommended for noisy environments (restaurants or social settings)
- Program 3** Recommended for quiet environments (1-on-1 conversations or watching TV)
- Program 4** T-Coil - Recommended for clearer sound when using a telephone or cell phone

##### Volume:

There are 5 volume levels on your Clarity Chat. A voice prompt will tell you when you have reached minimum or maximum volume.



### Troubleshooting your Clarity Chat

There are 2 probable issues that will cause you to believe your Clarity Chat isn't working properly: ear wax and low batteries.

#### Wax



See diagram 3 for cleaning instructions.

If wax buildup is in ear, please seek help from a doctor.

#### Low Battery



See diagram 1 for battery installation instructions.

#### Humidity


Other issues may be a result of humidity exposure. Humidity build up on the internal parts can cause the device to go dead or temporarily malfunction. If your device has been exposed to high humidity, remove it and open the battery door to let dry naturally.

#### Why is my Clarity Chat whistling?

Whistling without a tube on indicates your device is working properly. The main fix for whistling is to lower your volume 1-2 levels or reapply your Clarity Chat to your ear ensuring the dome is completely in the ear canal.


# DIAGRAM 1

## How to install the battery in your Clarity Chat:



# DIAGRAM 2


## How to wear your Clarity Chat:




# DIAGRAM 3

## How to change your tubing and domes:

### Tubing:



### Domes:



# DIAGRAM 4

## How to clean your Clarity Chat:

