vibrator is supplied as standard, then model, SS5V can be purchased as an option.

ONLY use the supplied Power Adaptor to operate this unit, the use of third party power adaptor

4) Attach the bed vibrator to the VIBRATOR INPUT (26) to the back of clock (if supplied). If no

OPERATING INSTRUCTION

1) Press and hold the TIME Button (22) and at the same time press the HOUR Button (18)

- b) In the mobile phone or tablet or mobile device, turn on the BT. Then search for "Sonic Blast Alarm Clock" and connect. Once connected, music can be played through the
- speaker of SB700. The volume of the music can be adjusted through the control of eh mobile device.
- To turn off the speaker, short press of the POWER button (13).
- 2) BT quick connection and disconnection

indicate it is searching for another BT device.

a) When the BT is connected, press VOL + (15) and VOL - (14) to connect (you will hear "connection" and BT light (5) light up). Press VOL + (15) and VOL - (14) again to disconnect (you will hear "disconnection" and BT light (5) flash).

Listening to Music

- There are three ways you can listen to music:
- 1) through BT 2) through SD card 3) through AUX input. To adjust the listening volume, press and hold the VOL - (14) button to decrease the
- volume. Press and hold the VOL + (15) button to increase the volume. 1) Music through BT
- a) Press and release the POWER button (13) to turn on the unit, press play on the mobile device to start playing music.
- 2) Music through Micro SD Card
- a) Plug the Micro SD card into SD socket (29). b) Press and release the POWER button (13) to turn on the unit. A short press of the MODE button (18) to select the "SD" mode. Then the music will start playing.
- <The below song selection applies only to music from BT and SD card>
- a) To select the next song, a short press of NEXT (15) button once; to select the previous song, a short press of the PREV (14) button once
- 3) Through AUX Input
- To connect the mobile device and the alarm clock, you need a 3.5mm to 3.5mm cable. a) Plug the one end of the 3.5mm cable into the AUX input (30) and the other end into
- b) Press and release the POWER button (13) to turn on this feature. A short press of the MODE button (18) to select the "AUX" mode, the AUX light (7) will light up.

- the mobile device (e.g MP3 player, mobile phone, ipod etc.)
- Then the music will start playing. Use mobile device to select songs.

- The default alarm duration time setting is 30 minutes; you may adjust the alarm duration between 1minute and 59 minutes as desired. To adjust the alarm duration: Press and hold the AL. SET button (23) for 6 seconds. The clock display will show the default alarm duration as '30'. Release the AL. SET button (23). 2) Press the HOUR (18) or MIN button (19) to adjust the alarm duration to the desired time,
- from 1 minute to 59 minutes. 3) It will return to the normal clock mode after 3 seconds of idle. The new alarm duration time

setting becomes the new default and remains in the memory until it is changed.

clock display will show the default snooze time '7'. Release the **SNOOZE Button (17)**.

becomes the new default and remains in the memory until it is changed.

NOTE: The dimmer function will not work, when the clock is in snooze mode

device to the USB charging port (28) on the back panel of this alarm clock.

1) To adjust the alarm volume, rotate VOLUME wheel (8).

Clock Display: Press the DIMMER button (17) once to change the display brightness.

There is a 5-level brightness control for the LED display.

2) To adjust the alarm Buzzer tone, turn the TONE wheel (11), the tone can be set from

This alarm clock is equipped with a dedicated USB charging port on the back of the clock

To charge your device, connect the USB charging cable that was supplied with your handled

a) When the clock is in standby (i.e. just display the clock, and nothing is on), press and

release the POWER button (13) to turn on this feature, the BT dot (5) will flash to

that allows you to recharge many brands of mobile phones, portable music players, and other

handheld devices that are supplied with USB charging cable. The output rating of the USB

2) Press the **MIN Button (19)** to adjust to the desired the snooze time, from 1 to 30 minutes.

3) It will return to the normal clock mode after 3 seconds of idle. The new snooze time setting

Adjusting The Snooze Time -

Dimmer Operation

Volume & Tone Control

800 - 1400Hz

charging port is 1A.

BT Speaker

1) Paring

USB Charging Function (DC5V, 1A)

Adjusting The Alarm Duration Time

or the MIN button (19) to advance clock display until you have desired time. The default snooze time is 7 minutes but you may adjust the snooze time between 1 and 30 2) The **lighted dot (1)** located in upper left corner of the clock display indicates the time is PM. minutes as desired. To adjust the snooze time: 1) Press and hold the SNÓOZE Button (17) for 3 seconds to enter the setting mode. The

Selecting Hours Format (12hr or 24 hr)

For AM time, this dot will not light up

Place vibrator between mattress and box spring

 Press and hold the HOUR Button (18) for 6 seconds to enter the hour selection mode, release the HOUR Button (18) when the display shows "12 H" or "24 H". 2) Short press the HOUR Button (18) to toggle between 12 or 24 hours.

Set the Alarm Time

is not recommended.

Set the Time

- Setting Alarm 1 1) Slide the ALARM SELECTOR (25) to AL1 position, the alarm indicator AL1 (2) will light up. 2) Press and hold the ALARM SET (23) button and at the same time press the HOUR Button (18) or the MIN button (19) to advance clock display until you have reach the desired
- alarm setting. 3) Make sure alarm time is correctly set for AM or PM (For 12 hour format).

Setting Alarm 2

Slide the ALARM SELECTOR (25) to AL2 position, the alarm indicator AL2 (3) will light up. Then repeat step 2 and 3 above to set alarm 2.

Alarm -Fnable Alarms

- 1) To enable Alarm 1 or Alarm 2, slide the ALARM SELECTOR (25) to their appropriate position. The alarm indicator AL1 (2) or AL2 (3) will light up accordingly.
- 2) To enable both alarms, slide the ALARM SELECTOR (25) to AL1+AL2 position, the alarm
- indicators AL1 (2) and AL2 (3) both will light up. 3) Slide the ALARM FUNCTION selector (10) to your desired wake up method, i.e. Vibrate only (VIB), Sound (ALM) or both.

Shutting Off the Alarm -

- Alarm can be shut off temporary for the day or permanently shut off. 1) When alarm goes off, press the POWER (13) button to turn off for the day and repeat at
- the same time the next day, if alarm settings have not changed. To shut the alarm off permanently, slide the ALARM FUNCTION (10) to the OFF position.
- NOTE: The Alarm Off function of Alarm 1 and Alarm 2 is independent of each other.

Snooze Operation -

- When the alarm turns on in the morning you may press the SNOOZE Button (17) to stop the alarm temporarily. The alarm will shut off for 7 minutes or for a predetermined time (see section under Adjusting The Snooze Time) and then turn on again. The snooze operation may be repeated indefinitely within 59 minutes.
- 2) The snooze function of Alarm 1 and Alarm 2 are independent of each other

To Select Wake up Alarm Sound

- To select the wake up alarm sound, slide the ALARM SOUND SELECTOR (9) to your desired
- a) AL. Play for Recording b) Sound 1 & 2 for internal pre-recorded sound
- c) BUZZER for buzzer sound

Using the Sleep Mode

- The Sleep mode lets you set the sleep timer for your music or internal sound. The sleep timer can be set from 15 to 90 minutes in 15 minutes interval.
- 1) A short press of the SLEEP button (19) to enter the SLEEP timer setting, Snooze dot (4) lighted up.
- 2) The default sleep timer setting is 90 minutes. A further short press of the SLEEP button
- (19) to change the timer time at interval of 15 minutes (90 / 60/ 45/ 30/ 15).
- 4) To cancel the Sleep Timer before it counts down to "0", a short press of the POWER button (13) to turn off the Sleep mode.
- 5) Once entered the sleep mode, a short press of the MODE (18) to select the below different play mode:
- a) internal pre-recorded sound b) BT (bluetooth) c) SD card d) AUX input 6) a) The internal pre-recorded playback volume can be adjusted by using the Alarm Volume Wheel (8)
- b) The BT, SD Card, AUX playback volume can be adjusted by using Vol (14) and Vol + (15)

SETUP AND POWER UP THE CLOCK Power Up this Unit

- 1) Unpack the unit, the DC power adapter and prepare 1 piece of CR2032 coin size battery (to be purchased separately). Remove all packing materials from the alarm clock.
- 2) Open the Battery cover (31). Insert the new battery into the compartment while observing the polarity and replace the battery compartment cover. (Backup battery is optional. Not installing this battery the unit will still work.)
- 3) Fully unwind the cord of the DC power adapter. Connect the male plug of the adapter to the DC5V Input (27) on the back of the clock. Insert the power adaptor into an AC outlet. We suggest that you use an outlet that is always 'live', not an outlet controlled by a wall switch

To prevent fire or shock hazard do not expose this appliance to rain or moisture

English

Model: SB700

• The ventilation should not be impeded by covering the ventilation openings with items such as, newspaper, table cloth, curtains, etc;

SAFETY INFORMATION

IMPORTANT - Please read these instructions carefully before use and retain for future reference.

- No naked flame sources, such as lighted candles should be placed on the apparatus.
- It shall not be exposed to dripping or splashing and that no objects filled with liquids shall be placed on the apparatus.

Only use the supplied or recommended Sonic Alert adaptor to connect the unit.

 Do not place the product in closed bookcases or racks without proper ventilation. . Where the mains/ power plug is used as the disconnect device, the disconnect device shall remain readily operable.

GENERAL DATA

4) Sleep Indicator5) BT Indicator

6) SD Indicator

12) Mic Input

7) AUX Indicator

8) Alarm Volume Wheel

9) Alarm Sound Selector

11) Buzzer Tone Wheel

10) Alarm Function Selector

210.00 mm


Warning:

UL listed Power Adaptor: AC100-240V~, 50/60Hz Power Consumption: Adaptor Output: DC5V ---- , 2200mA USB Charging: Max 1A SD Card: Max size supported 32GB (FAT32 format) Bluetooth: v4.2

Sonic Alert

Sonic Blast Alarm Clock with Projector

CONTROLS & INDICATORS



600000 ෂ්ත්ම මම ම -33

15) Next Track / Music Vol + 26) 5V vibrator input 27) Power Adaptor input DC5V 16) Alarm Strobe Light 28) USB Charging Port 17) Snooze / Display Dimmer 18) Mode / Hour 29) Micro SD Card slot 19) Sleep On / Min 30) AUX Input 20) Projector Lens 31) Backup Battery Compartment 21) Record /Record Playback 32) Projector Focus 33) Projection ON/Off / 22) Time Set 23) Alarm Set

Image Rotate

Clock Projection

1) Point the projection lens to the wall, press the Projection ON/OFF button (33) to turn on the projection function. Press the Projection Rotate button (33) to rotate the projected

2) Use the Focusing wheel (32) to focus the projection lens to make a clear projection image.

Voice Recording & Playback

- Voice Recording
 Press and hold the REC/PLAY button (21) for one second, then you will hear a beep to indicate that it has entered the recording mode.
- Speak to the MIC (12) to be recorded
- There will be two beeps to indicate it is the end of the recording. The recording can be up to a maximum of 20 seconds
- Voice Recording Playback

1) To playback the voice recording, a short press of REC/PLAY button (21).

Test Mode

The test mode is used to test the basic function of the clock

Hold the TIME (22) and SNOOZE (17) for 3 seconds to enter the test mode.
 Slide the FUNCTION SELECTOR (10) to test the different functions.

3) To exit press POWER button (13) or it will auto exit after 1 minute of inactivity.

FCC Statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference.

(2) This device must accept any interference received, including interference that may cause undesired operation

Note: The grantee is not responsible for any changes or modifications not expressly approved by the party responsible for compliance. Such modifications could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that

interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance of 20 cm between the radiator and your body

Normal function may be disturbed by strong electromagnetic interference. If this is the case, normal function can be resumed by simply resetting the unit following the procedure in the instructions manual. If function does not resume, try moving the unit to another location.

SONIC ALERT'S 1 YEAR LIMITED WARRANTY

The Sonic Blast Alarm Clock Model SB700 is warranted against manufacturing defects in materials and workmanship for one (1) year from the date of purchase.

Within this period Sonic Alert will repair or replace at our option the SB700 without charge for parts and labor. Simply contact our customer service department to obtain an RA

(return authorization number) and instructions on how to process your product for return. Call 1-888-864-2446 or e-mail to RMA@sonicalert.info.

You'll need a copy of your receipt to process the return authorization.

Sonic Alert Inc., Troy MI 48083 USA

Sonic Alert Sonic Alert, Troy, MI 48083 www.SonicAlert.com 1-888-864-2446

Made in China

093-SB700-final Rev02