






Written Programming Directions

STEP 1	STEP 2	STEP 3
<p>SETTING ALARM 1</p>  <p>Set Button</p> <p>Right Arrow Button</p> <ol style="list-style-type: none">1. Hold down the "SET" button & press the "RIGHT ARROW" button "1 TIME".2. The "HOUR" will start flashing. Use the "RIGHT ARROW" button & scroll to the desired hour and release.3. Press the "SET" button, the "MINUTES" will flash. Again, use the "RIGHT ARROW" button to scroll to your desired minute and release.4. Press the "SET" button, an "A OR P" will flash (A for AM & P for PM). Use the "RIGHT ARROW" button to select your desired time of day. <p>RELEASE & ALARM 1 IS SET</p>	<p>SETTING ALARM 2</p> <ol style="list-style-type: none">1. Hold down the "SET" button & press the "RIGHT ARROW" button "2 TIMES"2. The "HOUR" will start flashing. Use the "RIGHT ARROW" button & scroll to the desired hour and release.3. Press the "SET" button, the "MINUTES" will flash. Again, use the "RIGHT ARROW" button to scroll to your desired minute and release.4. Press the "SET" button, an "A OR P" will flash (A for AM & P for PM). Use the "RIGHT ARROW" button to select your desired time of day. <p>RELEASE & ALARM 2 IS SET</p> <hr/> <p>DISABLE ALARM 2</p> <p>This <u>MUST BE DONE</u>, when medication is taken ONCE DAILY.</p> <ol style="list-style-type: none">1. Hold down the "SET" button & press the "RIGHT ARROW" button "2 TIMES". The "HOUR" will start flashing. Use the "RIGHT ARROW" button to scroll past the #12 to the word "OFF" <p>RELEASE & ALARM 2 IS DISABLED</p>	<p>SETTING TIME & DAY</p> <ol style="list-style-type: none">1. Hold down the "SET" button & press the "RIGHT ARROW" button "3 TIMES"2. The "HOUR" will start flashing. Use the "RIGHT ARROW" button & scroll to the desired hour and release.3. Press the "SET" button, the "MINUTES" will flash. Again, use the "RIGHT ARROW" button to scroll to your desired minute and release.4. Press the "SET" button an "A OR P" will flash (A for AM & P for PM). Use the "RIGHT ARROW" button to select your desired time of day and release.5. Press the "SET" button and at the top of the screen 1 of the 7 days of the week will flash. (Compare to image below, inside the YELLOW rectangle). Use the "RIGHT ARROW" button to scroll to correct day.  <p>RELEASE & SET UP IS COMPLETE</p>

MORE INFORMATION

- After medication is taken, press the "NEXT ALERT" button to disarm alarm and activate the next alert.
- The Back row is either PM Pills OR WEEK 2.



- The Front row is AM Pills or WEEK 1.
- Escalating Alarm System
 - 100% volume for 2 min.
 - 150% volume for 2 min.
 - 200% volume for 1 min.
- Smart Sleep Mode
 - Hold down the "NEXT ALERT" button & "RIGHT ARROW" button.

WAKE UP: Hold down the "NEXT ALERT" button and "RIGHT ARROW" button. Wait 20 seconds, press "NEXT ALERT" button, now your alarms are back on.

MORE INFORMATION

MedqPillbox.com

602-488-4099

3919 W. Abraham Lane
Glendale, AZ 85308



The MED-Q is intended to serve as a backup reminder system to help remind users to take their medication. It is not intended to be the sole source for users to remember or properly take their medication. You agree that the limit of liability from any claim arising from the use, or malfunction of the MED-Q shall be limited to the purchase price.

Med-Q

INSTRUCTIONS



Please Visit
medqpillbox.com
For Video SET-UP
Instructions

**"Med-Q Remembers so YOU
DONT HAVE TO!"**

medqpillbox.com / 602-488-4099

Med-Q

Med-Q HOME

Set-Up Videos

Contact Us

Frequently Asked Questions

Buy Med-Q

Quick Contact

E: info@medqpillbox.com
Ph: 602 488-4099
Add: 3915 W Abraham Glendale,
AZ 85308

Med-Q