



## Presentation

The natural lubricant agent is strongly recommended for new hearing aid users.

Sweet almond oil makes it easier to insert the hearing aids. It also acts preventively to avoid irritation or discomfort concerns.

The natural lubricant agent is presented in a 0.50 fl oz (15mL) glass bottle, provided with an applicator brush.

## Properties

The natural lubricant meets the requirements of the Class I Medical Devices, according to the European Directive 93/42/EEC related to Medical Devices, while respecting environmental and safety requirements for health.

## Using advices

Apply the product on the earmold thanks to the brush for a better fitting and comfort in the ear.

## Precautionary statements

External use.  
Keep out of reach of children.  
Do not swallow.

## Product composition

- Sweet almond oil

